# Safety Data Sheet

Revision Date 2023-12-29

**Revision Number** 8

	1. Identification
Product identifier	
Product Name	50X Titanium Taq DNA Polymerase
Other means of identification	
Product Code	S1792
Synonyms	No information available
Recommended use of the chemic	al and restrictions on use
Identified uses	For research use only. Not for use in diagnostic procedures
Restrictions on use	No information available
Details of the supplier of the safe	ty data sheet
<u>Supplier</u> USA: Takara Bio USA, Inc. 2560 Orchard Parkway San Jose, CA 95131, USA Phone: 800.662.2566/888.251.6 Web: www.takarabio.com	3618
Emergency telephone number	
Emergency telephone	In case of emergency, call PERS (Professional Emergency Resource Services) 1-800-633-8253 (US) or 801-629-0667 (international).
	2. Hazard(s) identification
Product Classification Data	

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Acute toxicity - Inhalation (Dusts/Mists)

## Hazards not otherwise classified (HNOC)

Not applicable


## Label elements

Warning

## Hazard statements

This chemical is not considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200). Harmful if inhaled

Category 4



Avoid breathing dust/fume/gas/mist/vapors/spray Use only outdoors or in a well-ventilated area

## **Precautionary Statements - Response**

IF INHALED: Remove person to fresh air and keep comfortable for breathing Call a POISON CENTER or doctor if you feel unwell

#### Other information

No information available.

## 3. Composition/information on ingredients

#### Substance

Not applicable.

## Mixture

Chemical name	CAS No.	Weight-%	Trade secret
Glycerol	56-81-5	60 - 70	*

## 4. First-aid measures

## **Description of first aid measures**

General advice	Show this safety data sheet to the doctor in attendance.		
Inhalation	Remove to fresh air. If symptoms persist, call a physician. If breathing has stopped, give artificial respiration. Get medical attention immediately.		
Eye contact	Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.		
Skin contact	Wash skin with soap and water.		
Ingestion	Do NOT induce vomiting. Rinse mouth. Never give anything by mouth to an unconscious person. Get medical attention.		
Self-protection of the first aider	Ensure that medical personnel are aware of the material(s) involved, take precautions to protect themselves and prevent spread of contamination. Avoid breathing vapors or mists. Use personal protective equipment as required. See section 8 for more information.		

## Most important symptoms and effects, both acute and delayed

Symptoms Coughing and/ or wheezing. Difficulty in breathing.

## Indication of any immediate medical attention and special treatment needed

**Note to physicians** Treat symptomatically.

## 5. Fire-fighting measures

Suitable Extinguishing Media	Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.	
Large Fire	CAUTION: Use of water spray when fighting fire may be inefficient.	
Unsuitable extinguishing media	Do not scatter spilled material with high pressure water streams.	
Specific hazards arising from the chemical	No information available.	
Explosion Data Sensitivity to mechanical impac	t None.	
Sensitivity to static discharge	None.	
Special protective equipment and precautions for fire-fighters	Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear. Use personal protection equipment.	

## 6. Accidental release measures

#### Personal precautions, protective equipment and emergency procedures

Personal precautions	Ensure adequate ventilation. Avoid breathing vapors or mists. Use personal protective equipment as required.	
Other information	Refer to protective measures listed in Sections 7 and 8.	
Methods and material for containment and cleaning up		
Methods for containment	Prevent further leakage or spillage if safe to do so.	
Methods for cleaning up	Pick up and transfer to properly labeled containers.	

## 7. Handling and storage

#### Precautions for safe handling

Advice on safe handling Handle in accordance with good industrial hygiene and safety practice. Avoid breathing vapors or mists. Ensure adequate ventilation. In case of insufficient ventilation, wear suitable respiratory equipment. Do not eat, drink or smoke when using this product.

#### Conditions for safe storage, including any incompatibilities

**Storage Conditions** Keep containers tightly closed in a dry, cool and well-ventilated place. Keep out of the reach of children.

## 8. Exposure controls/personal protection

## Control parameters

## **Exposure Limits**

Chemical name	ACGIH TLV	OSHA PEL	NIOSH
Glycerol	-	TWA: 15 mg/m <sup>3</sup> mist, total	-
56-81-5		particulate	
		TWA: 5 mg/m <sup>3</sup> mist, respirable	
		fraction	
		(vacated) TWA: 10 mg/m <sup>3</sup> mist,	
		total particulate	
		(vacated) TWA: 5 mg/m <sup>3</sup> mist,	
		respirable fraction	

## Appropriate engineering controls

Engineering controls	Showers Eyewash stations Ventilation systems.
Individual protection measures, suc	ch as personal protective equipment
Eye/face protection	No special protective equipment required.
Skin and body protection	No special protective equipment required.
Respiratory protection	No protective equipment is needed under normal use conditions. If exposure limits are exceeded or irritation is experienced, ventilation and evacuation may be required.
General hygiene considerations	Do not eat, drink or smoke when using this product.

## 9. Physical and chemical properties

## Information on basic physical and chemical properties

Physical state	Liquid		
Appearance	Clear, colorless		
Color	No information available		
Odor	None		
Odor Threshold	No information available		
Property	Values	Remarks • Method	
pH	No data available	None known	
Melting point / freezing point	No data available	None known	
Boiling point/boiling range (°C)	No data available	None known	
Flash point	No data available	Open cup	
Evaporation Rate	No data available	None known	
Flammability (solid, gas)	No data available	None known	
Flammability Limit in Air		None known	
Upper flammability limit:	No data available		
Lower flammability limit:	No data available		
Vapor pressure	No data available	None known	
Vapor density	No data available	None known	

Relative density	No data available None kr			
Water solubility	No data available	None known		
Solubility in other solvents	No data available	None known		
Partition coefficient	No data available	None known		
Autoignition temperature	No data available	None known		
Decomposition temperature		None known		
Kinematic viscosity	No data available	None known		
Dynamic Viscosity	No data available	None known		
Other information				
Explosive properties	No information available			
Oxidizing properties	No information available			
Softening point	No information available			
Molecular weight	No information available			
VOC content	No information available			
Liquid Density	No information available			
Bulk Density	No information available			
10. Stability and reactivity				
		caotivity		
Reactivity	No information available.			

Reactivity	No mormation available.
Chemical stability	Stable under normal conditions.
Possibility of hazardous reactions	None under normal processing.
Conditions to Avoid	Excessive heat.
Incompatible materials	None known based on information supplied.
Hazardous decomposition products	None known based on information supplied.

## 11. Toxicological information

## Information on likely routes of exposure

## **Product Information**

Inhalation	Specific test data for the substance or mixture is not available. Harmful by inhalation. (based on components).		
Eye contact	Specific test data for the substance or mixture is not available.		
Skin contact	Specific test data for the substance or mixture is not available.		
Ingestion	Specific test data for the substance or mixture is not available.		
Symptoms related to the physical, chemical and toxicological characteristics			
Symptoms	Coughing and/ or wheezing.		
Acute toxicity			
Numerical measures of toxicity			
The following values are calculated based on chapter 3.1 of the GHS document			

#### Component Information

Chemical name	Oral LD50	Dermal LD50	Inhalation LC50		
Glycerol	= 12600 mg/kg (Rat)	> 10 g/kg (Rabbit)	> 2.75 mg/L (Rat)4 h		
56-81-5			i]		
Delayed and immediate effects as well as chronic effects from short and long-term exposure					
Skin corrosion/irritation	No information available.				
Serious eye damage/eye irritati	on No information available.				
Respiratory or skin sensitizatio	<b>n</b> No information available.				
Germ cell mutagenicity	No information available.				
Carcinogenicity	No information available.				
Reproductive toxicity	No information available.				
STOT - single exposure	No information available.				
STOT - repeated exposure	No information available.				
Aspiration hazard	No information available.				
Other adverse effects	No information available.				
Interactive effects	No information available.				

## 12. Ecological information

## Ecotoxicity

Chemical name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Glycerol 56-81-5	-	LC50: 51 - 57mL/L (96h, Oncorhynchus mykiss)	-	-

Persistence and degradability No in

No information available.

## **Bioaccumulation**

There is no data for this product.

## **Component Information**

Chemical name	Partition coefficient
Glycerol 56-81-5	-1.75

Other adverse effects

No information available.

13. Disposal considerations		
Waste treatment methods		
Waste from residues/unused products	Dispose of in accordance with local regulations. Dispose of waste in accordance with environmental legislation.	
Contaminated packaging	Do not reuse empty containers.	

## 14. Transport information

DOT

Not regulated

## 15. Regulatory information

#### International Inventories

## TSCA

\*Contact supplier for details. One or more substances in this product are either not listed on the US TSCA inventory, listed on the confidential US TSCA inventory or are otherwise exempted from inventory listing requirements

DSL/NDSL	
EINECS/ELINCS	
ENCS	
IECSC	
KECL	
PICCS	
AICS	

Legend:

**TSCA** - United States Toxic Substances Control Act Section 8(b) Inventory

DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List

**EINECS/ELINCS** - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances

**ENCS** - Japan Existing and New Chemical Substances

**IECSC** - China Inventory of Existing Chemical Substances

**KECL** - Korean Existing and Evaluated Chemical Substances

**PICCS** - Philippines Inventory of Chemicals and Chemical Substances

AICS - Australian Inventory of Chemical Substances

## US Federal Regulations

## SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product does not contain any chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372.

## SARA 311/312 Hazard Categories

Should this product meet EPCRA 311/312 Tier reporting criteria at 40 CFR 370, refer to Section 2 of this SDS for appropriate classifications.

#### CWA (Clean Water Act)

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42).

#### **CERCLA**

This material, as supplied, does not contain any substances regulated as hazardous substances under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302) or the Superfund Amendments and Reauthorization Act (SARA) (40 CFR 355). There may be specific reporting requirements at the local, regional, or state level pertaining to releases of this material.

#### US State Regulations

#### **California Proposition 65**

This product does not contain any Proposition 65 chemicals.

## U.S. State Right-to-Know Regulations

Chemical name	New Jersey	Massachusetts	Pennsylvania
Glycerol	X	Х	Х
56-81-5			

### U.S. EPA Label Information

## EPA Pesticide Registration Number Not applicable

16. Other information				
NFPA HMIS	Health hazards 2 Health hazards 2	Flammability 1 Flammability 1	Instability 0 Physical hazards 0	Special hazards - Personal protection X
Key or legend to abbreviations and acronyms used in the safety data sheetLegendSection 8: EXPOSURE CONTROLS/PERSONAL PROTECTIONTWATime weighted averageSTELCeilingMaximum limit value*Skin designation				
Agency for Toxic U.S. Environmer European Food EPA (Environmer Acute Exposure U.S. Environmer U.S. Environmer Food Research Hazardous Subs International Uni Japan GHS Clas Australia Nationa NIOSH (Nationa	stance Database iform Chemical Information Datab	ry (ATSDR) Database secticide, Fungicide, and uction Volume Chemicals ase (IUCLID) n and Assessment Scher and Health)	Rodenticide Act	

National Library of Medicine's PubMed database (NLM PUBMED) U.S. National Toxicology Program (NTP) New Zealand's Chemical Classification and Information Database (CCID) Organization for Economic Co-operation and Development Environment, Health, and Safety Publications Organization for Economic Co-operation and Development High Production Volume Chemicals Program Organization for Economic Co-operation and Development Screening Information Data Set World Health Organization

Revision Date	2023-12-29
Revision Note	No information available.

## **Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

**End of Safety Data Sheet**