

## **ICELL8® Loading Kit**

Catalog No.	Amount	Lot Number
640109 (Also sold as part of 640000, 640010, 640110,	1 pkg	Specified on product label.
640146, 640164, 640165 & 640166)		

## Description

The ICELL8 Loading Kit contains films and blotting paper designed to be used with ICELL8 chips. The Blotting Paper is used to blot chips after dispensing of cells into a chip by the ICELL8 MultiSample NanoDispenser (MSND). The Optical Imaging Film is used to seal chips prior to imaging of the nanowells with the ICELL8 Imaging System. The TE Sealing Film is used to seal the chip after reagents have been dispensed into cells selected for processing for reverse transcription.

## **Package Contents**

- 3 x Blotting Paper
- 1 x Optical Imaging Film
- 1 x TE Sealing Film

## **Storage Conditions**

• Store at room temperature

## Shelf Life

• 1 year from date of receipt under proper storage conditions.

## **Shipping Conditions**

• Room temperature

## **Product Documents**

Documents for our products are available for download at takarabio.com/manuals

## **Quality Control Data**

All elements of the ICELL8 Loading Kit have been visually inspected and found to be free of defects.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.



# **ICELL8® Loading Kit**

### CATALOG NO.

640109

### NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <u>http://www.takarabio.com</u>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

### **TRADEMARKS:**

### ©2017 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

<b>Takara Bio USA, Inc</b> . 1290 Terra Bella Avenue, Mountain View, CA 94043, USA U.S. Technical Support: techUS@takarabio.com					
United States/Canada	Asia Pacific	Europe	Japan	3/15/2019	
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999		