

Single-Cell Collection Kit

Catalog No.

Amount Each **Lot Number** Specified on product label.

640212 Also sold as part of 640188, 640257, 640282, 640286 & 640288

Description

The Single-Cell Collection Kit consists of a collection fixture, 2.0 ml tubes, and collection film. This kit can be used to spin down a processed Single-Cell 350v Chip to collect the contents (e.g., cDNA or library product) of the chip into a tube.

Package Contents

- 1 x Single-Cell Collection Fixture
- 2 x Collection Tube (2.0 ml)
- 1 x Collection Film

Storage Conditions

• Store at room temperature.

Shelf Life

• 1 year from date of receipt under proper storage conditions.

Shipping Conditions

• Room temperature

Product Documents

Documents for our products are available for download at <u>takarabio.com/manuals</u> The following documents apply to this product:

• ICELL8® cx Single-Cell System User Manual

Quality Control Data

All elements of the Single-Cell Collection Kit have been visually inspected and found to be free of defects.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.



Single-Cell Collection Kit

CATALOG NO.

640212

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <u>http://www.takarabio.com</u>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 436

This Product is protected by one or more patents from the family comprising: US9447925, People's Republic of China Patent: ZL201380058480.0, US10500585, DE602013083531.4, EP2916953, FR2916953, UK2916953, JP6559568, SE2916953 and any corresponding patents, divisionals, continuations, patent application and foreign filings sharing priority with the same family.

TRADEMARKS:

©2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.				
2560 Orchard Parkway, San Jose, CA 95131, USA				
U.S. Technical Support: technical_support@takarabio.com				0.15.1000.1
United States/Canada	Asia Pacific	Europe	Japan	8/5/2024
800.662.2566	+1.650.919.7300	+33.(0)1.3904.6880	+81.(0)77.565.6999	