

# Certificate of Analysis

---

## Single-Cell Extraction Fixture

**Catalog No.**

640217 (Not sold separately)  
Sold as a part of 640188 & 640282

**Amount**

Each

**Lot Number**

Specified on product label.

### Description

The Single-Cell Extraction Fixture is a device designed to hold the Single-Cell Collection Kit (Cat. No. 640212) during the process of extracting the components of a 5,184-nanowell Single-Cell chip into a single tube.

### Package Contents

- 1 Single-Cell Extraction Fixture

### Storage Conditions

- Store at room temperature.

### Shipping Conditions

- Room temperature

### Product Documents

Documents for our products are available for download at [takarabio.com/manuals](https://takarabio.com/manuals)

- Shasta™ Single Cell System User Manual
- ICELL8® cx Single-Cell System User Manual

### Quality Control Data

The Single-Cell Extraction Fixture has been visually inspected and found to be free of defects.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

---

**Takara Bio USA, Inc.**

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

United States/Canada  
800.662.2566  
(091924)

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.565.6999

## Single-Cell Extraction Fixture

### CATALOG NO.

640217

### NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

### STATEMENT 329

This Product is protected by one or more patents from the family consisting of: US8252581, CA2677833, US9132427, US9951381, US11643681, and any corresponding patents, divisionals, continuations, patent applications and foreign filings sharing common priority with the same family.

### TRADEMARKS:

#### ©2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

#### United States/Canada

800.662.2566

#### Asia Pacific

+1.650.919.7300

#### Europe

+33.(0)1.3904.6880

#### Japan

+81.(0)77.565.6999

9/19/2024