

Capturem™ Trypsin Activation Buffer

Catalog No.	Amount	Lot Number
635739	50 ml	Specified on product label.

Description

Capturem Trypsin Activation Buffer is used for activating Capturem Trypsin Miniprep Columns (Cat. No. 635740) and Capturem Trypsin 96-Well Plates (Mass Spectrometry Grade) (Cat. Nos. 635736, 635737). 50 ml of activation buffer is sufficient to activate two 96-well plates.

Package Contents

- 50 ml Capturem Trypsin Activation Buffer

Storage Conditions

- Room temperature

Shelf Life

- 1 year from date of receipt under proper storage conditions.

Shipping Conditions

- Room temperature

Product Documents

Documents for our products are available for download at takarabio.com/manuals

The following documents apply to this product:

- Capturem Trypsin 96-Well Plate (Mass Spectrometry Grade) Protocol-At-A-Glance

Quality Control Data

The Capturem Trypsin Activation Buffer has been tested for pH and conductivity.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada 800.662.2566 (053019)	Asia Pacific +1.650.919.7300	Europe +33.(0)1.3904.6880	Japan +81.(0)77.565.6999
--	---------------------------------	------------------------------	-----------------------------

Capturem™ Trypsin Activation Buffer

CATALOG NO.

635739

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 347

This product is protected by U.S. Patents 9,895,665, 10,195,569 and foreign counterparts and/or additional U.S. and foreign patents pending. For further license information, please contact a Takara Bio USA licensing representative by email at licensing@takarabio.com.

TRADEMARKS:

© 2019 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

1290 Terra Bella Avenue, Mountain View, CA 94043, USA

U.S. Technical Support: techUS@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

8/9/2019