

## Cellartis® CM Thawing Base

**Catalog No.**

Y10062

**Amount**

32 ml

**Lot Number**

Specified on product label

### Description

Cellartis CM Thawing Base is a medium used for thawing Cellartis Cardiomyocytes.

### Package Contents

- 32 ml Cellartis CM Thawing Base

### Storage Conditions

- Store at -20°C

### Shelf Life

- Specified on product label

### Shipping Conditions

- Dry ice (-70°C)

### Product Documents

Documents for Takara Bio Europe AB products are available for download at [www.takarabio.com/manuals](http://www.takarabio.com/manuals)

The following documents apply to this product:

- Cellartis Cardiomyocytes User Manual
- Cellartis Pure Cardiomyocytes User Manual

### Quality Control Data

#### Bacterial and Fungal Contamination Test

Cellartis CM Thawing Base has been tested and found to be free of bacterial and fungal contamination.

### Additional Notes

Use Cellartis CM Thawing Base with addition of 20% FBS and 10 µM Y-27632 to thaw Cellartis Cardiomyocytes and Cellartis Pure Cardiomyocytes.

### Notice to Purchaser

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

This product may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Europe AB.

If you require licenses for other use, please contact us by phone at +46.(0)31.758.0900.

Your use of this product is also subject to compliance with any applicable licensing requirements as detailed below, in our catalogues, on our website at <http://www.takarabio.com>, on the label or other documentation accompanying the goods. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

### TRADEMARKS

© 2018 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

---

**Takara Bio Europe AB**

A Takara Bio Company

Arvid Wallgrens backe 20, SE-413 46 Göteborg, Sweden,

Technical Support: [techEU@takarabio.com](mailto:techEU@takarabio.com)