Certificate of Analysis



Lenti-X™ Tet-Off® Advanced Vector Set

Catalog No. **Amount Lot Number** $10 \mu g$

632158 (Not sold separately) Sold as a part of 632163

Specified on product label.

Description

The Lenti-X Tet-Off Advanced Vector Set contains lentiviral expression vectors that are designed to establish an inducible expression system in mammalian cells. Lentivirus produced from the pLVX-Tet-Off Advanced Vector expresses a tetracycline-dependent transcriptional activator (tTA-Advanced) which specifically induces transcription of your gene of interest from the P_{Tight} inducible promoter in pLVX-Tight-Puro. In cells cotransduced by both lentiviral vectors, your cDNA is expressed to high levels in the absence of doxycycline or after its withdrawal. Using these vectors with the Lenti-X HTX Packaging System (Cat. No. 631247) will produce very high-titer supernatants from 293T packaging cells. The VSV-G pseudotyped lentiviruses so generated can be used to transduce virtually any mammalian cell type.

Package Contents

- 10 μg pLVX-Tet-Off Advanced Vector (500 ng/μl)
- 10 μg pLVX-Tight-Puro Vector (500 ng/μl)
- 10 μg pLVX-Tight-Puro-Luc Vector (500 ng/μl)

Storage Conditions

Store at -20°C.

Expiration Date

Specified on product label.

Storage Buffer

- 10 mM Tris-HCl (pH 8.0)
- 1 mM EDTA (pH 8.0)

Concentration

 $500 \text{ ng/}\mu\text{l}$

Shipping Conditions

Dry ice

Product Documents

Documents for our products are available for download at takarabio.com/manuals The following documents apply to this product:

- pLVX-Tet-Off Advanced Vector Information Packet
- pLVX-Tight-Puro Vector Information Packet
- pLVX-Tight-Puro-Luc Vector Information Packet

Quality Control Data

Plasmid Identity & Purity

• The identity of each plasmid was verified by sequencing critical elements, and by electrophoresis on a 0.8% agarose/EtBr gel after restriction enzyme digestion.:

Plasmid Name	Enzyme	Fragments (bp)
pLVX-Tet-Off Advanced	XhoI KpnI	11,578 9,718 & 1,860
pLVX-Tight-Puro	XhoI SalI	7,191 5,606 & 2,185
pLVX-Tight-Puro-Luc	XhoI SalI	9,424 7,239 & 2,185

• The purity of each plasmid was also checked by determining the A_{260}/A_{280} .

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

(121724) Page 2 of 2



Lenti-XTM Tet-Off® Advanced Vector Set

CATALOG NO.

632158

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at http://www.takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

STATEMENT 42

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. Patent # 8383364, # 9181556, European patents EP # 1954811, #2352833 and corresponding patent claims outside these regions which are proprietary to TET Systems GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions whose research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to any data and results stemming from the sponsored research, shall need a commercial license agreement from TET Systems in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Takara Bio USA, Inc. is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to TET Systems.

For license information, please contact: GSF/CEO
TET Systems GmbH & Co. KG,
Im Neuenheimer Feld 582
69120 Heidelberg
Germany

Tel: +49 6221 5880400 Fax: +49 6221 5880404 email: info@tetsystems.com

or use the electronic licensing request form via http://www.tetsystems.com/ip-licensing/licensing/for-profit-research

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA U.S. Technical Support: technical support@takarabio.com

United States/Canada Asia Pacific Europe Japan 12/17/2024

800.662.2566 +1.650.919.7300 +33.(0)1.3904.6880 +81.(0)77.565.6999

Notice to Purchaser



TRADEMARKS:

©2024 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

United States/Canada

2560 Orchard Parkway, San Jose, CA 95131, USA U.S. Technical Support: technical_support@takarabio.com

Asia Pacific

Europe Japan 800.662.2566 +1.650.919.7300 +81.(0)77.565.6999 +33.(0)1.3904.6880