

## Magnetic Beads Immunoprecipitation Buffer Set

**Catalog No.**  
635696

**Amount**  
50 rxns

**Lot Number**  
Specified on product label.

### Description

The Magnetic Beads Immunoprecipitation Buffer Set is a set of lysis, wash, and elution buffers optimized for purifying proteins immunoprecipitated with Anti-DYKDDDDK Magnetic Beads or Anti-Myc Magnetic Beads.

### Package Contents

- 110 ml Lysis Buffer
- 2 x 100 ml Wash Buffer
- 2 ml IP Elution Buffer

### Storage Conditions

- Store at 4°C.

### Expiration Date

- Specified on product label.

### Shipping Conditions

- Blue ice

### Product Documents

Documents for our products are available for download at [takarabio.com/manuals](http://takarabio.com/manuals)

The following documents apply to this product:

- Anti-Myc Immunoprecipitation Protocol-At-A-Glance
- Magnetic DYKDDDDK Immunoprecipitation Kit Protocol-At-A-Glance

### Quality Control Data

The buffers were tested for pH and conductivity.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

United States/Canada  
800.662.2566  
(062823)

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.565.6999

## Magnetic Beads Immunoprecipitation Buffer Set

### CATALOG NO.

635696

### NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

### TRADEMARKS:

©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

---

#### Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: [technical\\_support@takarabio.com](mailto:technical_support@takarabio.com)

#### United States/Canada

800.662.2566

#### Asia Pacific

+1.650.919.7300

#### Europe

+33.(0)1.3904.6880

#### Japan

+81.(0)77.565.6999

6/28/2023