


pTRE-Tight-BI-ZsGreen1 Vector Information

PT3875-5

Catalog No. 631067

Visit our website for more details! click here...



		<u>SmaI</u>		<u>NheI</u>
	<u>KpnI</u>		<u>BamHI</u>	
331	GCTCGGTACC	CGGGGATCCT	CTAGTCAGCT	GACGCGTGCT
	CGAGCCATGG	GCCCCTAGGA	GATCAGTCGA	CTGCGCACGA
		<u>NotI</u>	<u>HindIII</u>	<u>Sall</u>
	<u>NheI</u>	<u>EagI</u>	<u>Clal</u>	<u>AccI</u>
371	AGCGCGGCCG	CATCGATAAG	CTTGTCGACG	
	TCGCGCCGGC	GTAGCTATTC	GAACAGCTGC	

pTRE-Tight-BI-ZsGreen1 Vector Map and Multiple Cloning Site (MCS). All sites shown are unique.

Description

pTRE-Tight-BI-ZsGreen1 is a bidirectional TRE-Tight plasmid that can be used to inducibly express a green fluorescent protein (ZsGreen1) along with a gene of interest with our Tet-On® and Tet-Off® Gene Expression Systems and Cell Lines (1, 2). The Tet Expression Systems and Cell Lines provide researchers ready access to the tetracycline-regulated expression systems described by Gossen & Bujard (3; Tet-Off) and Gossen *et al.* (4; Tet-On).

pTRE-Tight-BI-ZsGreen1 contains a modified Tet response element (TRE_{mod}), which consists of seven direct repeats of a 36 bp sequence that contains the 19 bp tet operator sequence (*tetO*) (5; pTREtight). The two mini CMV promoters (P_{minCMV}), which lack the enhancer that is part of the complete CMV promoter, flank the TRE_{mod}. pTRE-Tight-BI-ZsGreen1 encodes a variant of wild-type *Zoanthus sp.* green fluorescent protein (ZsGreen1) (excitation maximum = 493 nm; emission maximum = 505 nm). pTRE-Tight-BI-ZsGreen1 contains a multiple cloning site (MCS) downstream of the bidirectional, Tet-responsive $P_{Tight-BI}$ promoter. Both ZsGreen1 and genes inserted into the MCS will be responsive to the tTA and rTA regulatory proteins in the

(PR0Z3756; published 6 December 2010)



Clontech

United States/Canada
800.662.2566

Asia Pacific
+1.650.919.7300

Europe
+33.(0)1.3904.6880

Japan
+81.(0)77.543.6116

Clontech Laboratories, Inc.
A Takara Bio Company
1290 Terra Bella Ave.
Mountain View, CA 94043
Technical Support (US)
E-mail: tech@clontech.com
www.clontech.com

Tet-Off and Tet-On systems, respectively. Note that the cloned insert must have an initiation codon (ATG). In some cases, addition of a Kozak consensus sequence (6) may improve expression levels; however, many cDNAs have been efficiently expressed in Tet systems without the addition of a Kozak sequence. pTRE-Tight-BI-ZsGreen1 should be cotransfected with a Linear Hygromycin or Puromycin Marker (Cat. Nos. 631625 and 631626, respectively; not included) to permit selection of stable transfectants (7). pTRE-Tight-BI-ZsGreen1 was derived from pTRE, (originally described as pUHD10-3 (8)) and pTREtight (5).

The pTRE-Tight-BI-Luc Control Vector, packaged with the pTRE-Tight-BI-ZsGreen1 Vector, lacks the ZsGreen1 gene, but contains a 1,649 bp firefly luciferase gene inserted into the MCS. This vector can be used as a reporter of induction efficiency using standard luciferase detection reagents.

Location of features

- $P_{\text{Tight-BI}}$ (bidirectional, Tet-responsive promoter):
 - TRE_{mod} (modified Tet-response element): 3–252
 - $P_{\text{minCMV-1}}$ (minimal CMV promoter 1): 258–317
 - $P_{\text{minCMV-2}}$ (minimal CMV promoter 2): 3440–3509 (complementary)
- MCS (multiple cloning site): 335–399
- SV40 polyA signal: 417–617
- pUC origin of replication: 780–1379
- Amp^r (ampicillin resistance gene; β -lactamase): 1540–2536 (complementary)
- SV40 polyA signal: 2538–2728 (complementary)
- ZsGreen1 gene: 2734–3429 (complementary)

Propagation in *E. coli*

- Suitable host strains: DH5 α^{TM} and other general purpose strains.
- Selectable marker: plasmid confers resistance to ampicillin (100 $\mu\text{g/ml}$) in *E. coli* hosts.
- *E. coli* replication origin: pUC

Excitation and emission maxima of ZsGreen1

- Excitation maximum = 493 nm
- Emission maximum = 505 nm

References

1. HT1080 Cell Line & pTRE2 Vector (January 1999) *Clontechniques* **XIV**(1):23.
2. Matz, M.V. *et al* (1999) *Nature Biotech.* **17**(10):969-973
3. Gossen, M. & Bujard, H. (1992) *Proc. Natl. Acad. Sci USA* **89**(12):5547–5551.
4. Gossen, M. *et al.* (1995) *Science* **268**(5218):1766–1769.
5. pTREtight vectors (April 2003) *Clontechniques* **XVIII** (2):10–11.
6. Kozak, M. (1987) *Nucleic Acids Res.* **15**(20):8125–8148.
7. Linear Selection Markers (April 2003) *Clontechniques* **XVIII**(2): 11.
8. Resnitzky, D. *et al.* (1994) *Mol. Cell. Biol.* **14**:1669–1679.

Notice to Purchaser

Clontech products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without written approval of Clontech Laboratories, Inc.

DH5 α [™] is a trademark of Life Technologies Corporation.

The RCFP's (including DsRedExpress and DsRedExpress2) are covered by one or more of the following U.S. Patent Nos. 7,166,444; 7,157,565; 7,217,789; 7,338,784; 7,338,783; 7,537,915 6,969,597, 7,150,979 and 7,442,522.

Not-For-Profit Entities: Orders may be placed in the normal manner by contacting your local representative or Clontech Customer Service at 650.919.7300. At its discretion, Clontech grants Not-For-Profit Entities a non-exclusive, personal, limited license to use this product for non-commercial life science research use only. Such license specifically excludes the right to sell or otherwise transfer this product, its components or derivatives thereof to third parties. No modifications to the protein coding sequence may be made without express written permission from Clontech. Any other use of this product requires a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com.

For-Profit Entities wishing to use this product are required to obtain a license from Clontech. For license information, please contact a licensing representative by phone at 650.919.7320 or by e-mail at licensing@clontech.com or click here for more information.

Use of the Tetracycline controllable expression systems (the "Tet Technology") is covered by a series of patents including U.S. Patent Nos. 5,464,758 and 5,814,618, which are proprietary to TET Systems GmbH & Co. KG. Academic research institutions are granted an automatic license with the purchase of this product to use the Tet Technology only for internal, academic research purposes, which license specifically excludes the right to sell, or otherwise transfer, the Tet Technology or its component parts to third parties. Notwithstanding the above, academic and not-for profit research institutions whose research using the Tet Technology is sponsored by for profit organizations, which shall receive ownership to all data and results stemming from the sponsored research, shall need a commercial license agreement from TET Systems in order to use the Tet Technology. In accepting this license, all users acknowledge that the Tet Technology is experimental in nature. TET Systems GmbH & Co. KG makes no warranties, express or implied or of any kind, and hereby disclaims any warranties, representations, or guarantees of any kind as to the Tet Technology, patents, or products. All others are invited to request a license from TET Systems GmbH & Co. KG prior to purchasing these reagents or using them for any purpose. Clontech is required by its licensing agreement to submit a report of all purchasers of the Tet-controllable expression system to TET Systems. For license information, please contact: GSF/CEO, TET Systems GmbH & Co. KG, Im Neuenheimer Feld 582, 69120 Heidelberg, Germany Tel: +4962215880400, Fax: +4962215880404 eMail: info@tetsystems.com or use the electronic licensing request form via http://www.tetsystems.com/main_inquiry.htm

Clontech, the Clontech logo and all other trademarks are the property of Clontech Laboratories, Inc., unless noted otherwise. Clontech is a Takara Bio Company. ©2010 Clontech Laboratories, Inc.