

PrimePath™ Saliva Collection Kit

Catalog No.

638343

Also sold as a part of 638340

Amount

200 preps

Lot Number

Specified on product label.

Description

The PrimePath Saliva Collection Kit contains collection sticks and a viral inactivation buffer that is used both for storage and processing samples. At one end of each collection stick, there are channels for collecting liquid such as saliva. When mixed with the collected saliva sample, the buffer will inactivate any microbes that are present. This process makes the saliva compatible with the RT-qPCR step in the PrimePath Saliva Direct SARS-CoV-2 Detection Kit.

Package Contents

- 100 packs NextSwab Universal Sterile Applicator (2 sticks/pack)
- 20 ml PrimePath Viral Inactivation Buffer

Storage Conditions

- Store at Room temperature

Expiration Date

- Specified on product label.

Shipping Conditions

- Room temperature

Product Documents

Documents for our products are available for download at [takarabio.com/manuals](https://www.takarabio.com/manuals)

The following documents apply to this product:

- PrimePath Direct Saliva SARS-CoV-2 Detection Kit User Manual

Quality Control Data

Physical tests:

The PrimePath Viral Inactivation Buffer was tested for pH, conductivity, density, and nuclease contamination.

It is certified that this product meets the above specifications, as reviewed and approved by the Quality Department.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA

U.S. Technical Support: technical_support@takarabio.comUnited States/Canada
800.662.2566
(042523)Asia Pacific
+1.650.919.7300Europe
+33.(0)1.3904.6880Japan
+81.(0)77.565.6999

PrimePath™ Saliva Collection Kit

CATALOG NO.

638343

NOTICE TO PURCHASER:

Our products are to be used for **Research Use Only**. They may not be used for any other purpose, including, but not limited to, use in humans, therapeutic or diagnostic use, or commercial use of any kind. Our products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without our prior written approval.

Your use of this product is also subject to compliance with the licensing requirements, listed below if applicable, and described on the product's web page at <http://www.takarabio.com>. It is your responsibility to review, understand and adhere to any restrictions imposed by these statements.

TRADEMARKS:

©2023 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.

Takara Bio USA, Inc.

2560 Orchard Parkway, San Jose, CA 95131, USA
U.S. Technical Support: technical_support@takarabio.com

United States/Canada

800.662.2566

Asia Pacific

+1.650.919.7300

Europe

+33.(0)1.3904.6880

Japan

+81.(0)77.565.6999

4/25/2023