

## Cellartis® iPSC Single-Cell Cloning DEF-CS™ Additives

**Catalog No.**

Y30019 (Not sold separately)  
Sold as a part of Y30021

**Amount**

1 kit

**Lot Number**

Specified on product label

### Description

The Cellartis iPSC Single-Cell Cloning DEF-CS Additives product is sold as a part of the Cellartis iPSC Single-Cell Cloning DEF-CS Culture Media Kit, a defined culture system for efficient generation of human pluripotent stem cell clones from seeded single cells. In addition to the additives, the complete system includes basal medium and coating substrate.

### Package Contents

- 2 x 750 µl DEF-CS GF-1
- 500 µl DEF-CS GF-2
- 500 µl DEF-CS GF-3

### Storage Conditions

- Store at –20°C.

### Shelf Life

- Specified on product label

### Shipping Conditions

- Dry ice

### Product Documents

Documents for Takara Bio Europe AB products are available for download at [www.takarabio.com/manuals](http://www.takarabio.com/manuals)  
The following documents apply to this product:

- Cellartis iPSC Single-Cell Cloning DEF-CS Culture Media Kit User Manual

### Additional Notes

The Cellartis iPSC Single-Cell Cloning DEF-CS Additives product is designed to be used with Cellartis DEF-CS 500 Basal Medium (Cat. No. Y30011; not sold separately) and Cellartis iPSC Single-Cell Cloning DEF-CS COAT-1 (Cat. No. Y30018; not sold separately). The basal medium, coating substrate, and additives are components of the Cellartis iPSC Single-Cell Cloning DEF-CS Culture Media Kit (Cat. No. Y30021).

### Quality Control Data

#### Bacterial and Fungal Contamination Test

DEF-CS GF-1, DEF-CS GF-2, and DEF-CS GF-3 have been tested and found to be free of bacterial and fungal contamination.

---

**Takara Bio Europe AB**

A Takara Bio Company

Arvid Wallgrens backe 20, SE-413 46 Göteborg, Sweden

Technical Support: [techEU@takarabio.com](mailto:techEU@takarabio.com)

United States/Canada  
800.662.2566  
(022322)

Asia Pacific  
+1.650.919.7300

Europe  
+33.(0)1.3904.6880

Japan  
+81.(0)77.565.6999

# Certificate of Analysis

Cat. No. Y30019

Cellartis iPSC Single-Cell Cloning DEF-CS Additives (Not sold separately)

Sold as a part of Y30021

---

## Notice to Purchaser

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

This product may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Europe AB.

If you require licenses for other use, please contact us by phone at +46 (0)31.758.0900.

Your use of this product is also subject to compliance with any applicable licensing requirements as detailed below, in our catalogues, on our website at <http://www.takarabio.com>, on the label or other documentation accompanying the goods. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

## TRADEMARKS

© 2019 Takara Bio Inc. All Rights Reserved.

All trademarks are the property of Takara Bio Inc. or its affiliate(s) in the U.S. and/or other countries or their respective owners. Certain trademarks may not be registered in all jurisdictions.